

## LAB ALERT

## Date: June 28, 2016

## HIV and H. Pylori Testing Updates

Dear Regional Pathology Services Clients,

We are very excited to announce upcoming changes in our serology testing that will tremendously improve test turn-around times. Effective July 1<sup>st</sup>, 2016, the Microbiology Laboratory will be switching to the BioPlex 2200 platform which utilizes 5<sup>th</sup> generation HIV testing that detects **HIV-1 antibody**, **HIV-2 antibody**, and **HIV-1 p24 antigen** each of which will be reported as individual results.

The HIV-1 p24 antigen is useful in the detection of HIV infection during the window period of 3 months following acute infection.

All positive HIV antibody and antigen results will be confirmed. This new HIV Antigen/Antibody Panel (HIVPN) will replace HIV1 and 2 Antibody ELISA Screen (HIV12). Specimen requirements have not changed.

Code CPT Code	Result Code	Result Name	Old Reference Range	New reference Range
N 87389, 87390, 86701, 86702	HIVAA HIV1AB HIV1AG HIV2AB	HIV AG-AB HIV1 Ab HIV1 p24 Ag HIV2 Ab	Non-Reactive Reactive	No Change
	Code N 87389, 87390, 86701,	Code Code   N 87389, 87390, 86701, HIVAA HIV1AB HIV1AG	CodeCodeN87389, 87390, 86701,HIVAA HIV1AB HIV1AGHIV AG-AB HIV1 Ab HIV1 Ab HIV1 p24 Ag	CodeCodeRangeN87389, 87390, 86701,HIVAA HIV1AB HIV1AGHIVAG-AB HIV1Ab HIV1 p24 AgNon-Reactive Reactive

Interface Clients: Please deactivate HIV12 order code and add the HIVPN code

*II. Helicobacter pylori* IgM, IgG and IgA testing discontinued effective July 1<sup>st</sup>.

Due to lack of sensitivity, the clinical laboratory is no longer offering *Helicobacter pylori* serology testing. The *H. pylori* stool antigen test or urease breath test is recommended for diagnosis of this infection. *Interface clients: Please deactivate order code HPYAB*.

For any questions about these changes please contact client services at 402-559-6420 and ask to speak with one of our client coordinators.

24/7 Client Services 1-800-334-0459